CHONDROGRID

PRODUKTINFORMATION

KORTIKOIDE

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HYDROLYSIERTES KOLLAGEN

ANALGETIKA

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
ChondroGrid - Studien

Article Intra-Articular Injection of Hydrolyzed Collagen to Treat Symptoms of Knee Osteoarthritis. A Functional In Vitro Investigation and a Pilot Retrospective Clinical Study

Paola De Luca ^{1,†}, Alessandra Colombini ^{1,†}, Giulia Carimati ², Michelangelo Beggio ³, Laura de Girolamo ^{1,*} and Piero Volpi ²

¹ IRCCS Istituto Ortopedico Galeazzi, Via R. Galeazzi 4, 20161 Milano, Italy
² Istituto Clinico Humanitas, Via Alessandro Manzoni, 56, 20089 Rozzano Milano, Italy
³ Policlinico San Marco, Via Francesco Zanotto, 40, 30173 Mestre, Italy

* Correspondence: laura.degirolamo@grupposandonato.it; Tel.: +39-0266214069

† These authors contributed equally to this work.

Received: 30 May 2019; Accepted: 3 July 2019; Published: 4 July 2019

Abstract: Among all joints affected, knee osteoarthritis has a prevalence of about 10% in men and 13% in women over 60 years old. Knee osteoarthritis has high economic and social costs and may have a devastating impact on patient quality of life. Treatment of symptomatic knee Osteoarthritis may involve oral or topical administration of non-steroidal anti-inflammatory drugs or intra-articular injection of corticosteroids. Recently, a novel injectable collagen formulation (ChondroGrid) consisting of bovine hydrolyzed <3 kDa type I collagen has been developed and is currently available on the market as an injectable medical device. The primary objective of this study was to investigate the in vitro and in vivo effects of ChondroGrid in treating knee osteoarthritis symptoms to assess its safety and performance. Viability and proliferation of ChondroGrid-exposed human chondrocytes derived from five donors were assessed through the Alamar Blue/CyQuant assays. Their expression of MMP1/MMP3 and TIMP1/TIMP3 was then assessed through RT-PCR and that of TGFβ1, IGF-I, and VEGF using ELISA assays. Shape and ECM deposition were assessed using the Bern score after a 28-day ChondroGrid exposure, and collagen deposition was assessed using immunostaining. Records of 20 patients affected by Kellgren Lawrence grade 1 to 4 knee osteoarthritis who received three 4 mg/2 mL ChondroGrid injections 2 weeks apart were then retrospectively assessed to compare VAS, Lequesne, and WOMAC scores collected before and 15, 45, and 225 days after the first injection. ChondroGrid had no effects on the markers under consideration, but induced type-II and inhibited type-I collagen deposition; the Bern score was higher when cells were cultured with ChondroGrid. Patients experienced a 44% Lequesne score and a 55% VAS at moving score reduction. All other scores decreased >70%. ChondroGrid may prompt chondrocytes to produce hyaline cartilage, prevent fibrous tissue formation, and be a safe and effective adjuvant to treat symptomatic knee osteoarthritis.

Keywords: knee; osteoarthritis; hydrolyzed collagen; intra-articular injection; non-pharmacological therapy

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International Orthopaedics
<https://doi.org/10.1007/s00264-020-04616-8>

ORIGINAL PAPER

Effectiveness of a novel hydrolyzed collagen formulation in treating patients with symptomatic knee osteoarthritis: a multicentric retrospective clinical study

Piero Volpi ¹ · Raul Zini ² · Franz Erschbaumer ³ · Michelangelo Beggio ⁴ · Alberto Busilacchi ² · Giulia Carimati ¹

Received: 30 January 2020 / Accepted: 11 May 2020
 © SICOT aisbl 2020

Abstract

Purpose: Knee osteoarthritis (OA) is a musculoskeletal disorder that may have a heavy impact on the patients' quality of life. Intra-articular collagen injection may be a safe adjuvant. Recently, CHondroGrid (CG), a hydrolyzed (< 3 kDa) bovine collagen injectable formulation, has been placed on the market. The aim of this study was to investigate the safety and performance profile of CG.

Methods: Patients affected by Kellgren Lawrence grade 1 to 4 knee OA and BMI < 30 were treated by administering three CG injections of 2 ml (4 mg) each (at 15 days and 45 days from the first one, respectively) and were followed up for six months after the last administration. Clinical records were retrospectively assessed to compare VAS, Lequesne and WOMAC total, pain, stiffness, and physical function scores collected at baseline and 15, 45, and 225 days after the first injection.

Results: At the last follow-up, 70 patients (37 men and 33 women, aged 57.1 ± 14.5 years) treated with CG showed a 50% reduction in their median Lequesne score, a 50% reduction in their VAS score at rest and moving, and a ≥ 50% reduction for all other scores under consideration.

Conclusions: CG may be a safe and effective adjuvant in the treatment of symptomatic knee OA.

Keywords: CHondroGrid · Intra-articular injection · Non-pharmacological therapy · Knee osteoarthritis

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Eurofins Product Testing Italy S.r.l. certifies that the manufacturer has submitted to evaluation the technical design documentation relating to the medical device listed above and that it has been successfully evaluated according to the procedure of Annex II point 4 of the Directive 93-42-EEC and updated.

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BIOTECK S.p.A.